Fit Test C & B

Instructions for use

Product description:

Fit Test C & B is an addition-curing silicone for checking the precision of the fit of crowns, bridges, inlays/onlays, secondary parts and attachments, as well as for temporary blocking-out during final cementing of dentures. For the dentist **Fit Test C & B** is a diagnostic help for the optimum fitting of the finished restorative piece. With its blue colouring it allows to clearly visualise faulty areas on metal and full ceramic restorative pieces.

Fit Test C & B convinces with its simple handling and particularly high tear resistance.

Indications:

Control of precise fit, marginal tightness and marginal gap of
 - full or partial crowns
 - bridges
 - inlays/onlays
 - cast post abutments/post crowns
 - metal bases
 - implant-supported crowns and bridges
 - attachments and other secondary parts of dentures

Application:

Clean and dry restoration, preparation and abutment as usual. Due to technical reasons slight differences in the filling level between the syringe chambers can occur prior to the first use. Therefore remove syringe cap and squeeze the material out until the material is delivered evenly from both chambers. Match the nooks on the syringe and the mixing tip and fix the mixing tip by a 90° turn clockwise. The material is mixed automatically in the correct ratio when dispensing. All QM syringes come with a detachable piston. The piston can simply be pulled from the QM syringe and re-placed, without causing any leakage of material. The plugs in the syringe chambers will stay in place. This prevents any cross-contamination of the material in the syringe.

Fit Test C & B is applied directly from the mixing tip or, in case of inaccessible areas, with an instrument to the cleaned restorative piece or cavity. Cover all surfaces with **Fit Test C & B** and avoid trapping air. Place the restorative piece in situ and fix (do **not** let the patient bite down!).

The working time from the beginning of the mixing is at least 1 min; the intraoral setting time is 2 min maximum – monitor setting by checking excess material. After setting carefully remove the restorative piece with the attached silicone film. Small one-surface inlays can also be removed with sticky wax.

Evaluation:

Even silicone layers show a good precision of fit. In case of uneven silicone layers, the restorative piece or the preparation has to be corrected using standard techniques in translucent contact areas.

Information, precautionary measures:

Fit Test C & B contains platinum catalyst.

Fit Test C & B should not be used in case of known hypersensitivities (allergies) to this ingredient.

Do not use **Fit Test C & B** on existing permanently soft, silicone-based base materials (e.g. liners) since a chemical bond to these materials may occur.

The following materials are able to delay or prevent the setting reaction of **Fit Test C & B**: condensation-curing silicones, polysulfide impression materials, eugenol-containing materials, incompletely cured synthetic materials, as well as latex gloves. When handling **Fit Test C & B** these materials must not be used. Avoid contamination (e.g. of textiles) with non-cured silicone since it may be difficult to remove. The above mentioned working times apply to a temperature of approx. 23 °C and a relative humidity of approx. 50 %. Lower temperatures prolong these times, higher temperatures shorten them.

The warranted product characteristics can only be achieved when using the original Sulzer Mixpac mixing tips as approved by **VOCO**. All the mixing tips for use with our respective products are available from **VOCO**.

Storage:

Store at temperatures of 4 °C to 23 °C. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA2095 E1 1217 99 © by VOCO
VOCO GmbH