


CPS PRIME START: A NEW BEGINNING

A holder for every person

Try the two holders and see which one suits you the best.




Simply click out of the packaging

The interdental brushes are extremely well protected, and can be clicked straight out of the packaging. It's incredibly simple and hygienic.

Extract without opening

Extracting a brush is very simple: position the holder, click on the brush, and pull out. Done.



The packaging is Recyclable

Both the brushes and the holder can be disposed of with your household rubbish. The packaging is now recyclable, however: one part can be recycled as PET, the other as cardboard. A little service to the environment.


The click system

Whether you prefer the UHS 409 or the UHS 470 holder, both options work for any brush thanks to the click system, meaning you only have to buy replacement brushes. This also means less waste.

Swiss Made

Just as the prime refill, this products will also be made in Switzerland.

 **MADE IN SWITZERLAND**



CPS prime brushes are highly popular: durably made, with resilient brushes for a wide cleaning radius and a slim wire for narrow interdental spaces. The refill packaging is environmentally friendly and offers excellent protection from dirt and dust. Available in five CPS prime sizes, all compatible with all holder types, such as the UHS 409 and UHS 470 holders. Time to change your interdental brush? Simply click it on and get started.



Instructions for clicking on the brush



UHS holder 409 click system



UHS holder 470 click system



▶ Instructional video
scan code or:
curaprox.com/refill



CURAPROX is the oral health brand of CURADEN AG, based in Kriens, Switzerland. Since 1972, CURAPROX has been working together with experts from research, teaching and practice to develop gentle and effective tools for oral hygiene.

www.curaprox.com

CURAPROX

✚ SWISS PREMIUM ORAL CARE