Imprint[™] 4 Bite

VPS Bite Registration Material



As fast and accurate as it can be

Imprint™ 4 Bite VPS Material offers you:

- ► Good flow behaviour For high accuracy
- ► Fast intra-oral setting
 For a faster procedure
- ► High final hardness and low flexibility

For easy cutting and trimming and accurate stone model positioning

- ➤ Fresh material colour and scannability For good readability and easy access to the digital workflow
- ► Cartridge delivery

 For simple and convenient mixing and dispensing

Developed for highly accurate bite registrations.

With its great flow properties and a high enough viscosity, Imprint 4 Bite material does not flow deeply into interdental areas and delivers highly accurate bite registrations. A high final hardness and low flexibility after setting allow for easy cutting and trimming as well as accurate positioning of stone models.

Saves time with a fast intra-oral setting.

Imprint 4 Bite offers a pleasant working time and a short setting time of only 60 seconds, keeping procedures moving quickly while enhancing patient comfort. And, it helps to reduce inaccuracies due to jaw movements during setting.

	Mixing and application	Setting time from start of mixing*
Imprint™ 4 Bite VPS Bite Registration Material	0:20	1:00

^{*}The stated times are valid at 23°C/74°F. Times depend on ambient temperature.



Indications:

Universal VPS bite registration material - ideally suited for occlusal recording including optical registration of occlusal data for CAD/CAM systems.



Get full control of your impression work.

3M ESPE impression solutions are designed to add speed, convenience and precision to the entire procedure.



Imprint™ 4 Bite VPS Material is part of the successful family of Imprint 4 VPS impression materials. Try also the versatile range of Imprint™ 4 VPS materials for your everyday precision impressions as well as Imprint 4 Preliminary VPS Material.

Ordering Information



Imprint™ 4 Bite Refill (71529)

Item No. **Product Information**

71529 Imprint™ 4 Bite Refill

2 cartridges - 50ml each, 5 Garant™ Mixing Tips green,

5 Garant™ Intra-oral Tips Bite

www.3MESPE.com



3M Deutschland GmbH Location Seefeld 3M ESPE · ESPE Platz 82229 Seefeld · Germany info3mespe@mmm.com www.3MESPE.com

3M, ESPE, Garant and Imprint are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.

© 2014, 3M. All rights reserved.

XX (0X.2014)